

kardiomed[®] 530

boosted user manual

General



Read the instruction manual with all the security advices and warnings **before the initial use** of the units in order to secure a safe and conventional use. Keep this document throughout the life time of the machine to be able to refer to the instructions. In case of resale, pass it on with the unit.

Overview

All kardiomed® 530 units are available as a boosted version. Those units are equipped with special elements to make the units even more robust. Besides that the units offer more comfort.

Unit	Max. user weight	Features
motion boosted cycle 800 fit/ med	200 kg	different elements to boost
motion boosted cross 800 fit/ med	200 kg	different elements to boost
motion boosted body 800 fit/ med	200 kg	different elements to boost
motion boosted relax 800 fit/ med	250 kg	Turnable seat, higher back rest, arm rest and different elements to boost
motion boosted sprint 800 fit/ med	250 kg	different elements to boost

Maintenance

The required safety level can only be ensured by checking the unit regularly for damage and attrition. Damaged parts need to get replaced right away and the unit needs to get decommissioned until its repair.

Check the unit visually before every use, but at least daily for damaged parts, loose connections and atypical noises.

Depending on the usage but the latest after six months the unit has to get inspected intensely. All connections and weld seams need to get checked and the unit needs to get cleaned.

Clean the devices with mild purifier and a wet cloth to get off aggressive sweat. Dry the device afterwards. The use of spray bottles is not allowed; it has shown that the unit never dries at spots that can't be reached easily, so that rust may occur. Pay attention that no liquids run into the monitor as it could destroy parts of the electronic.

If disinfection is required only use acryl-des® disinfection cloths (see <http://www.schuelke.com/>).

The monitor is showing a sign for low battery voltage when the functioning of the machine cannot be guaranteed without electricity from the brake anymore. Exchange the 3AA batteries by loosening the monitor from the frame. Take out the two screws on the back of the monitor and open it. Please pay attention that insert the new batteries the right way around, note the signing. Dispose the batteries according to the provisions. There is a return obligation; the customer can choose to give the batteries to the local disposer or to the distributor free of charge.

Only use original spare parts, therefore please contact us.

Warning: Retorque of cranks and pedals

After the first 10 hours in usage the cranks and pedals need to get retorquened.

If they do not get retightened, damage of the axis, pedals or cranks is possible and the warranty will not cover the resulting damage.

Warranty conditions

The following is based on a legal warranty.

The warranty is terminated as soon as the product is serviced or repaired by an unauthorised person. If a warranty case occurs call the proxomed® Service Hotline +49 (6023) 9168 77 immediately.

proxomed® will initiate service, but reserves the right to decide about the mode of service.

The following procedures are possible:

1. The service is conducted on-site by our service employees.
2. The required spare part(s) are shipped.
3. An exchange unit is arranged and sent to the user.

The defective parts should be sent to EMOTION FITNESS within 48 hours. Otherwise charges for spare parts that are delivered will be charged.

If causes for repairs lie outside of the warranty conditions, proxomed® reserves the right to require payment for all costs of the repair.

Some consumable parts are excluded from the warranty.

Manufacturer according to MPG
emotion fitness GmbH & Co KG
Trippstadter Str. 68
67691 Hochspeyer
Telefon: +49 (0)6305-71499-0
Telefax: +49 (0)6305-71499-111
E-Mail: info@emotion-fitness.de

Exklusiv sales of the proxomed® kardiomed 530 units

proxomed® Medizintechnik GmbH
Daimlerstraße 6
63755 Alzenau

Telefon : +49 (0)6023 9168-0
Telefax : +49 (0)6023 916868

Internet : www.proxomed.de
E-mail : info@proxomed.com

Head Office Germany:

proxomed® Medizintechnik GmbH

Daimlerstraße 6

D-63755 Alzenau

Tel.: +49 6023 9168-0

Fax: +49 6023 9168-68

www.proxomed.de

info@proxomed.de

Branch Office Switzerland:

proxomed® Medizintechnik

Seestrasse 161

CH-8266 Steckborn

Tel.: +41 52762 1300

Fax: +41 52762 1470

www.proxomed.ch

info@proxomed.ch

Edits:

proxomed® reserves the right to change products if the actions lead in our view to improvements in quality and functions. All pictures are solely attached due to printing issues; We don't assume liability for typing errors. Errors excepted.

Service manual kardiomed® 530, V 2.0 / 19.08.2014
Number of article:



proxomed® is certified through
TÜV SÜD Product Service GmbH.